

HemOxy12+ OXYGEN CONCENTRATORS

(A Make in India Initiative "Atma Nirbhar")

CE & ISO: 13485: 2016 CERTIFIED

Acute breathing problem is one of the major issues of patient. The patient's oxygen saturation level becomes a critical parameter in saving the life. Therefore, oxygen concentrators can be the best support system for the patients and in case of lack of proper transportation facility for oxygen cylinders this can really be useful in hospitals and even for home use.

Why **HemOxy** by HEMIND® are better Concentrators...

- Engineered to suit Indian climate
- Robust yet portable
- Consistency of purity of Oxygen in continuous & long use duty cycle
- Minimum maintenance
- Spares available at ease as made in India.
- Two channels 6LPM x 2 channels = 12 LPM with effective 93 % (+/- 3%) uninterrupted pure oxygen flow. 90% + purity at 14LPM, 82%+ at 16 LPM

Features of the PSA design:

1. Compact design (portable and point of use facility)
2. User-friendly Interface with oil- free compressor
3. Consistent purity output for longer duration.
4. Flow rate range: (Two channel.)

Model: **HemOxy12**: 12 LPM (Continuous duty cycle)



USAGE

Hospitals ~ Home ~ Offices ~ Underground Bunkers
Hotel Rooms ~ Oxygen Bars
Beauty Parlours ~ Spas

We are a Class I India Manufacturer with more than 80 % localized components.

Manufactured under License from ISRO. Tested CIC (R&D centre – DST approved)

For Purchase inquiries please contact:



HEMIND®

HEMANT INDUSTRIES

AN ISO 9001: 2015 COMPANY

269, GIDC ESTATE, MAKARPURA ROAD, BARODA 390010

CONTACT : +91 9825014516. Email : info@hemind.com. Website : www.hemind.com